

---

Product Manual

# Nuclear/Cytosolic Fractionation Kit


Catalog Number

AKR-171

20 preps

**FOR RESEARCH USE ONLY**  
**Not for use in diagnostic procedures**

---



**CELL BIOLABS, INC.**  
*Creating Solutions for Life Science Research*

## **Introduction**

Cell Biolabs' Nuclear/Cytosolic Fractionation Kit provides a simple and fast tool to isolate nuclear extract from the cytoplasmic fraction of mammalian cells. The procedure has been optimized to provide extraction, with high protein recovery and low cross-contamination, in less than 2 hours. The extracted protein fractions are functional and suitable for downstream assays such as DNA footprinting, RNA splicing, gel shift assays (EMSA), reporter assays, enzyme activity assays, and Western blotting. Each kit provides sufficient quantities to perform 20 preps (up to  $5 \times 10^6$  cells each).

## **Related Products**

1. AKR-103: PhosphoBlocker™ Blocking Reagent (1L)
2. AKR-172: Nuclear/Cytosolic Fractionation Kit (100 Preps)

## **Kit Components (shipped on blue ice packs)**

1. Cytosol Extraction Buffer, Hypotonic (10X) (Part No 217101): One 2 mL vial.
2. Cell Lysis Reagent (Part No. 217102): One 1 mL vial of 10% Igepal CA-630 in 1X Cytosol Extraction Buffer (CEB).
3. Nuclear Extraction Buffer (Part No. 217103): One 2 mL vial.
4. Dithiothreitol (1000X) (Part No. 217104): One 20 µL vial of 1 M DTT.
5. Protease Inhibitor Cocktail (100X) (Part No. 217105): One 200 µL vial containing AEBSF, Aprotinin, Bestatin, E64, Leupeptin, and Pepstatin A in DMSO.

## **Materials Not Supplied**

1. PBS
2. Microcentrifuge tubes
3. Microcentrifuge

## **Storage**

Upon receiving, aliquot and store Dithiothreitol and Protease Inhibitor Cocktail at -20°C and avoid multiple freeze/thaw cycles. Store all other components at 4°C.

## **Preparation of Reagents**

- 1X Cytosol Extraction Buffer (CEB): Dilute the 10X Cytosol Extraction Buffer to 1X with deionized water. Stir to homogeneity.
- Dithiothreitol: Immediately before use dilute the Dithiothreitol 1:1000 with 1X Cytosol or Nuclear Extraction Buffer. Stir to homogeneity. Do not store diluted solutions.
- Protease Inhibitor Cocktail: Immediately before use dilute the Protease Inhibitor Cocktail 1:100 with 1X Cytosol or Nuclear Extraction Buffer. Stir to homogeneity. Do not store diluted solutions.

## **Preparation of Samples**

### **I. Adherent Cells**

1. Culture cells to approximately 80-90% confluence.
2. Aspirate the culture media and wash twice with PBS.
3. Detach the cells from the plates in PBS by scraping with a cell scraper.
4. Collect the solution into an appropriate conical centrifuge tube.
5. Centrifuge for 5 minutes (600 x g).
6. Discard the supernatant and immediately proceed to the Assay Protocol.

### **II. Suspension Cells**

1. Collect the cells into an appropriate conical centrifuge tube.
2. Centrifuge for 5 minutes (600 x g).
3. Remove and discard the supernatant.
4. Wash the cells twice with PBS.
5. Centrifuge for 5 minutes at (600 x g).
6. Discard the supernatant and immediately proceed to the Assay Protocol.

## **Assay Protocol**

*Important Note: Perform the below steps at 2-8°C. All buffers, centrifuge rotors, and equipment should be maintained at 2-8°C. Before use, Dithiothreitol and Protease Inhibitor Cocktail should be diluted according to the Preparation of Reagents section above.*

### **I. Cytosol Fractionation Protocol**

1. Collect cells (up to  $5 \times 10^6$ ) by centrifugation for 5 minutes at 4°C (600 x g).
2. Wash the cells once with ice cold PBS.
3. Remove and discard the supernatant.
4. Gently resuspend the cell pellet with 500  $\mu$ L of ice cold, 1X Cytosol Extraction Buffer (containing DTT/Protease Inhibitors) by pipetting up and down.
5. Transfer the suspension into a prechilled microcentrifuge tube.
6. Incubate on ice for 10 minutes.
7. Add 25  $\mu$ L of Cell Lysis Reagent and vortex for 10 seconds at the highest setting.
8. Centrifuge for 10 minutes at 4°C (800 x g).
9. Carefully transfer the supernatant (cytoplasmic fraction) to a clean, chilled microcentrifuge tube. The cytoplasmic fraction can be stored at -80°C for future use.

*Note: Make sure not to disturb/remove the nuclei pellet.*

10. Gently resuspend the pellet with 500  $\mu$ L of ice cold, 1X Cytosol Extraction Buffer (containing DTT/Protease Inhibitors) by pipetting up and down.

*Note: This wash step is included to reduce cross-contamination between fractions.*

11. Add 25  $\mu$ L of Cell Lysis Reagent and vortex for 10 seconds at the highest setting.
12. Centrifuge for 10 minutes at 4°C (800 x g).
13. Carefully aspirate the supernatant and discard of this wash.

## **II. Nuclear Protein Extraction Protocol**

1. Gently resuspend the nuclear pellet with 100  $\mu$ L of ice cold, 1X Nuclear Extraction Buffer (containing DTT/Protease Inhibitors) by pipetting up and down.
2. Maintain on ice for 30 minutes, vortexing for 10 seconds at the highest setting in 10 minute intervals.
3. Centrifuge for 30 minutes at 4°C (14000 x g).
4. Carefully transfer the supernatant (nuclear protein extract) to a clean, chilled microcentrifuge tube. The extract can be stored at -80°C for future use.


*Note: The nuclear extract typically yields protein concentrations of > 1 mg/mL. If greater concentrations are desired, resuspend the nuclear pellet in a smaller volume in step 1 above (minimum of 25  $\mu$ L).*

## **III. Other Considerations**

- For determining the protein content of extracts, samples must be diluted 1:2 before running in the Bradford Protein Assay. Buffer only controls must be performed concurrently. DTT in the buffers is not compatible with the BCA Protein Assay.
- Nuclear Extraction Buffer is a high salt buffer, containing 420 mM NaCl. If salt removal is necessary, dialysis or a desalting column may be used.


## Example of Results

The following figure demonstrates typical results seen with Cell Biolabs' Nuclear/Cytosolic Fractionation Kit. One should use the data below for reference only.



**Figure 1: HEK293 Cell Fractionation.** Cytosolic and nuclear protein extracts were isolated from Human Embryonic Kidney 293 cells according to the Assay Protocol. Whole cell (W), cytosol (C), and nuclear (N) fractions were immunoblotted with Anti- $\alpha$ -Tubulin (left) or Anti-Lamin A/C (right) at 1  $\mu$ g/mL.

*Note: Anti- $\alpha$ -Tubulin (Calbiochem CP06) and Anti-Lamin A/C (Sigma SAB4200236) are both mouse monoclonals. Tubulin and Lamin are known to be cytosolic and nuclear specific proteins, respectively.*



**Figure 2: HEK293 Trypan Blue Staining.** Human Embryonic Kidney 293 cells were stained with Trypan Blue at various steps during the fractionation protocol, demonstrating complete lysis and high nuclei recovery.

## Recent Product Citations

1. Wang, X. et al. (2026). A ESRP1/circPHGDH/miR-149/RAP1B positive feedback loop promotes the malignant behaviors and glycolysis of prostate cancer cell. *Exp Mol Med.* **58**(2):622-635. doi: 10.1038/s12276-026-01646-x.

2. Kim, Y.L. et al. (2025). Tyr-Tyr-Glu tripeptide for regeneration of articular cartilage through chondrogenic differentiation of mesenchymal stem cells. *Biomater Res*. doi: 10.34133/bmr.0272.
3. Honda, N. et al. (2025). Implications of Mutant SOD1 on RNA Processing and Interferon Responses in Amyotrophic Lateral Sclerosis: Omics Data Analysis. *Cureus*. **17**(3):e81045. doi: 10.7759/cureus.81045.
4. Izkovich, B. et al. (2025). Virally mediated expression of a biologically active peptide to restrain the nuclear functions of ERK1/2 attenuates learning extinction but not acquisition. *Mol Brain*. **18**(1):19. doi: 10.1186/s13041-025-01190-1.
5. Fujimura, T. et al. (2024). Targeting ErbB and tankyrase1/2 prevent the emergence of drug-tolerant persister cells in ALK-positive lung cancer. *NPJ Precis Oncol*. **8**(1):264. doi: 10.1038/s41698-024-00757-w.
6. Tripathi, B.K. et al. (2024). The pro-oncogenic noncanonical activity of a RAS•GTP:RanGAP1 complex facilitates nuclear protein export. *Nat Cancer*. **5**(12):1902-1918. doi: 10.1038/s43018-024-00847-5.
7. Yamada, S. et al. (2024). SARS-CoV-2 causes dysfunction in human iPSC-derived brain microvascular endothelial cells potentially by modulating the Wnt signaling pathway. *Fluids Barriers CNS*. **21**(1):32. doi: 10.1186/s12987-024-00533-9.
8. Kim, H.J. et al. (2024). Low-Molecular-Weight Hydrolysate from Black Goat Extract Had Antioxidative and Anti-Inflammatory Effects in Macrophage Cells via Inhibition of MAPKs and NF- $\kappa$ B Pathways. *J Food Biochem*. doi: 10.1155/2024/7155015.
9. Moriya, K. et al. (2023). Human RELA dominant-negative mutations underlie type I interferonopathy with autoinflammation and autoimmunity. *J Exp Med*. **220**(9): e20212276. doi: 10.1084/jem.20212276.
10. Onishi, K. et al. (2023). Substrate stiffness induces nuclear localization of myosin regulatory light chain to suppress apoptosis. *FEBS Lett*. **597**(5):643-656. doi: 10.1002/1873-3468.14592.
11. Takagi-Kimura, M. et al. (2022). BAP1 depletion in human B-lymphoblast cells affects the production of innate immune cytokines and chemokines. *Genes Cells*. **27**(12):731-740. doi: 10.1111/gtc.12988.
12. Mohri, S. et al. (2022). Integration of bioassay and non-target metabolite analysis of tomato reveals that  $\beta$ -carotene and lycopene activate the adiponectin signaling pathway, including AMPK phosphorylation. *PLoS One*. **17**(7):e0267248. doi: 10.1371/journal.pone.0267248.
13. Li, Y. et al. (2022). LncRNA SNHG5 promotes the proliferation and cancer stem cell-like properties of HCC by regulating UPF1 and Wnt-signaling pathway. *Cancer Gene Ther*. doi: 10.1038/s41417-022-00456-3.
14. Mikawa, M. et al. (2022). Herpud1 suppress angiotensin II induced hypertrophy in cardiomyocytes. *Biochem Biophys Rep*. **30**:101248. doi: 10.1016/j.bbrep.2022.101248.
15. Nakazawa, N. et al. (2022). Cytoplasmic localization of connexin 26 suppresses transition of  $\beta$ -catenin into the nucleus in intestinal- and mix-type gastric cancer. *Ann Gastroenterol Surg*. doi: 10.1002/ags3.12552.
16. Xu, J. et al. (2022). The circular RNA circ\_0030018/miR-136/migration and invasion enhancer 1 (MIEN1) axis promotes the progression of polycystic ovary syndrome. *Bioengineered*. **13**(3):5999-6011. doi: 10.1080/21655979.2022.2041796.
17. Fu, Y. et al. (2022). Long non-coding RNA HCG22 inhibits the proliferation, invasion and migration of oral squamous cell carcinoma cells by downregulating miR-425-5p expression. *Exp Ther Med*. **23**:246. doi: 10.3892/etm.2022.11171.

18. Tripathi, B.K. et al. (2021). Inhibition of cytoplasmic EZH2 induces antitumor activity through stabilization of the DLC1 tumor suppressor protein. *Nat Commun.* **12**(1):6941. doi: 10.1038/s41467-021-26993-3.
19. Huang, D. & Li, C. (2021). circ-ACACA promotes proliferation, invasion, migration and glycolysis of cervical cancer cells by targeting the miR-582-5p/ERO1A signaling axis. *Oncol Lett.* **22**(5):795. doi: 10.3892/ol.2021.13056.
20. Takeuchi, K. et al. (2021). Colchicine protects against cartilage degeneration by inhibiting MMP13 expression via PLC- $\gamma$ 1 phosphorylation. *Osteoarthritis Cartilage.* doi: 10.1016/j.joca.2021.08.001.
21. Xu, Z. et al. (2021). Lnc-HZ01 with m6A RNA methylation inhibits human trophoblast cell proliferation and induces miscarriage by up-regulating BPDE-activated lnc-HZ01/MXD1 positive feedback loop. *Sci Total Environ.* doi: 10.1016/j.scitotenv.2021.145950.
22. Shu, J. et al. (2021). Fertility-enhancing potential of ethanol extract of *Cuscuta chinensis* seeds in a rat model of unilateral cryptorchidism. *Trop J Pharm Res.* **20**(5):995-1002. doi: 10.4314/tjpr.v20i5.16.
23. Geng, C. & Ou, S. (2021). Higenamine Induces Glioma Cell Death by Modulating Nuclear Factor- $\kappa$ B Nuclear Translocation, Phosphoinositide-3-Kinase/Protein Kinase B Signaling and Caspase Cascade. *Curr Top Nutraceutical Res.* **19**(3):317-325. doi: 10.37290/ctnr2641-452X.19:317-325.
24. Li, J. et al. (2021). LncRNA NOP14-AS1 Promotes Tongue Squamous Cell Carcinoma Progression by Targeting MicroRNA-665/HMGB3 Axis. *Cancer Manag Res.* **13**:2821-2834. doi: 10.2147/CMAR.S293322.
25. Bhatt, A.B. et al. (2021). Diverse and converging roles of ERK1/2 and ERK5 pathways on mesenchymal to epithelial transition in breast cancer. *Transl Oncol.* **14**(6):101046. doi: 10.1016/j.tranon.2021.101046.
26. Ushimaru, S. et al. (2020). Roles of Layilin in Regulation of Low-Density Lipoprotein Receptor in Malignant Glioma Cells. *J. St. Marianna Univ.* **11**:53-59.
27. Kumar, A. et al (2020). Actin R256 Mono-methylation Is a Conserved Post-translational Modification Involved in Transcription. *Cell Rep.* **32**(13):108172. doi: 10.1016/j.celrep.2020.108172.
28. Hwang, S.G. et al. (2020). Cold atmospheric plasma prevents wrinkle formation via an anti-aging process. *Plasma Med.* doi: 10.1615/PlasmaMed.2020034810.
29. Sadek, J. et al. (2020). Modulation of virus-induced NF- $\kappa$ B signaling by NEMO coiled coil mimics. *Nat Commun.* **11**(1):1786. doi: 10.1038/s41467-020-15576-3.
30. Xiao, G. et al. (2019). Bacoside A attenuates nephrotoxicity and acute kidney injury in male albino rats induced by cisplatin. *Int. J. Pharmacol.* **15**:257-264. doi: 10.3923/ijp.2019.257.264.

## **Warranty**

These products are warranted to perform as described in their labeling and in Cell Biolabs literature when used in accordance with their instructions. THERE ARE NO WARRANTIES THAT EXTEND BEYOND THIS EXPRESSED WARRANTY AND CELL BIOLABS DISCLAIMS ANY IMPLIED WARRANTY OF MERCHANTABILITY OR WARRANTY OF FITNESS FOR PARTICULAR PURPOSE. CELL BIOLABS' sole obligation and purchaser's exclusive remedy for breach of this warranty shall be, at the option of CELL BIOLABS, to repair or replace the products. In no event shall CELL BIOLABS be liable for any proximate, incidental or consequential damages in connection with the products.

## **Contact Information**

Cell Biolabs, Inc.  
5628 Copley Drive  
San Diego, CA 92111  
Worldwide: +1 858-271-6500  
USA Toll-Free: 1-888-CBL-0505  
E-mail: [tech@cellbiolabs.com](mailto:tech@cellbiolabs.com)  
[www.cellbiolabs.com](http://www.cellbiolabs.com)

©2012-2026: Cell Biolabs, Inc. - All rights reserved. No part of these works may be reproduced in any form without permissions in writing.