

## 293 LTV Cell Line

---

**CATALOG NUMBER:** LTV-100

**STORAGE:** Liquid nitrogen

*Note: For best results begin culture of cells immediately upon receipt. If this is not possible, store at -80°C until first culture. Store subsequent cultured cells long term in liquid nitrogen.*

**QUANTITY & CONCENTRATION:** 1 mL, >1 x 10<sup>6</sup> cells/mL in 90% complete medium, 10% DMSO

### **Background**

The 293LTV Cell Line is a permanent line established from primary embryonic human kidney transformed with human adenovirus type 5 DNA. The genes encoded by the E1 region of adenovirus (E1a and E1b) are expressed in these cells and participate in transactivation of viral promoters, allowing these cells to produce high levels of protein.

293LTV also stably expresses the SV40 large T antigen and Neomycin resistance gene, through cloning and multiple rounds of testing for viral yield, 293LTV is specifically selected for high level of lentiviral production. It offers several advantages over the regular 293T cells:

- High lentiviral yield
- Firm attachment to culture plate and fast growing
- Ideal as a host when making lentivirus by cotransfection

### **Quality Control**

This cryovial contains at least 1.0 × 10<sup>6</sup> 293LTV cells as determined by morphology, trypan-blue dye exclusion, and viable cell count. The 293LTV cells are tested free of microbial contamination.

### **Medium**

1. Culture Medium: D-MEM (high glucose), 10% fetal bovine serum (FBS), 0.1 mM MEM Non-Essential Amino Acids (NEAA), 2 mM L-glutamine, 1% Pen-Strep (optional)
2. Freeze Medium: 90% complete medium, 10% DMSO

### **Methods**

#### **I. Establishing 293LTV Cultures from Frozen Cells**

1. Place 10 mL of complete DMEM growth medium in a 50-mL conical tube. Thaw the frozen cryovial of cells within 1–2 minutes by gentle agitation in a 37°C water bath. Decontaminate the cryovial by wiping the surface of the vial with 70% (v/v) ethanol.
2. Transfer the thawed cell suspension to the conical tube containing 10 ml of growth medium.
3. Collect the cells by centrifugation at 1000 rpm for 5 minutes at room temperature. Remove the growth medium by aspiration.
4. Resuspend the cells in the conical tube in 15 mL of fresh growth medium by gently pipetting up and down.

5. Transfer the 15 mL of cell suspension to a T-75 tissue culture flask. Place the cells in a 37°C incubator at 5% CO<sub>2</sub>.
6. Monitor cell density daily. Cells should be passaged when the culture reaches 95% confluence.

### **Recent Product Citations**

1. Long, F. et al. (2026). MEDAG functions as an A-kinase-anchoring protein in adipocytes. *Mol Cell*. **86**(5):937-953.e9. doi: 10.1016/j.molcel.2026.02.001.
2. Whitford, A.L. et al. (2025). Interferon dependent immune memory during HSV-1 neuronal latency via increased H3K9me3 and restriction by ATRX. *Nat Commun*. doi: 10.1038/s41467-025-66477-2.
3. Wu, C. et al. (2025). MTCH2 modulates CPT1 activity to regulate lipid metabolism of adipocytes. *Nat Commun*. **16**(1):8831. doi: 10.1038/s41467-025-63880-7.
4. Liu, Y.C. et al. (2025). APP SUMOylation prevents BACE1 cleavage of APP and increases BACE1 degradation to promote the nonamyloidogenic pathway. *Mol Med*. **31**(1):301. doi: 10.1186/s10020-025-01354-8.
5. Asamori, T. et al. (2025). Molecular mimicry-driven autoimmunity in chronic rhinosinusitis with nasal polyps. *J Allergy Clin Immunol*. doi: 10.1016/j.jaci.2025.02.014.
6. Lai, C.Y. et al. (2024). CtBP1 is essential for epigenetic silencing of  $\mu$ -opioid receptor genes in the dorsal root ganglion in spinal nerve ligation-induced neuropathic pain. *Neurotherapeutics*. **22**(1):e00493. doi: 10.1016/j.neurot.2024.e00493.
7. Wang, T. et al. (2024). Single-nucleus transcriptomics identifies separate classes of UCP1 and futile cycle adipocytes. *Cell Metab*. **36**(9):2130-2145.e7. doi: 10.1016/j.cmet.2024.07.005.
8. Safarzadeh Kozani, P. et al. (2024). Humanization of the antigen-recognition domain does not impinge on the antigen-binding, cytokine secretion, and antitumor reactivity of humanized nanobody-based CD19-redirected CAR-T cells. *J Transl Med*. **22**(1):679. doi: 10.1186/s12967-024-05461-8.
9. Cheng, K.M. et al. (2024). Novel role of NCoR1 in impairing spatial memory through the mediation of a novel interacting protein DEC2. *Cell Mol Life Sci*. **81**(1):273. doi: 10.1007/s00018-024-05321-0.
10. Hilliard W. et al. (2023). A compendium of stable hotspots in the CHO genome. *Biotechnol Bioeng*. doi: 10.1002/bit.28390.
11. Kuhse, J. et al. (2023). Loss of Extrasynaptic Inhibitory Glycine Receptors in the Hippocampus of an AD Mouse Model Is Restored by Treatment with Artesunate. *Int J Mol Sci*. **24**(5):4623. doi: 10.3390/ijms24054623.
12. Proaño-Pérez, E. et al. (2023). MITF Downregulation Induces Death in Human Mast Cell Leukemia Cells and Impairs IgE-Dependent Degranulation. *Int J Mol Sci*. **24**(4):3515. doi: 10.3390/ijms24043515.
13. Vullhorst, D. et al. (2023). ER-PM Junctions on GABAergic Interneurons Are Organized by Neuregulin 2/VAP Interactions and Regulated by NMDA Receptors. *Int J Mol Sci*. **24**(3):2908. doi: 10.3390/ijms24032908.
14. Liu, B. et al. (2022). Functional genomics screening identifies aspartyl-tRNA synthetase as a novel prognostic marker and a therapeutic target for gastric cancers. *J Pathol*. doi: 10.1002/path.5980.
15. Moy, R.H. et al. (2022). Functional genetic screen identifies ITPR3/calcium/RELB axis as a driver of colorectal cancer metastatic liver colonization. *Dev Cell*. **57**(9):1146-1159.e7. doi: 10.1016/j.devcel.2022.04.010.

16. Fish, L. et al. (2021). A prometastatic splicing program regulated by SNRPA1 interactions with structured RNA elements. *Science*. **372**(6543):eabc7531. doi: 10.1126/science.abc7531.
17. Balazova, L. et al. (2021). GPR180 is a component of TGF $\beta$  signalling that promotes thermogenic adipocyte function and mediates the metabolic effects of the adipocyte-secreted factor CTHRC1. *Nat Commun*. **12**(1):7144. doi: 10.1038/s41467-021-27442-x.
18. Navinés-Ferrer, A. et al. (2021). MYO1F Regulates IgE and MRGPRX2-Dependent Mast Cell Exocytosis. *J Immunol*. **206**(10):2277-2289. doi: 10.4049/jimmunol.2001211.
19. Wade, M. et al. (2021). In vivo generation of collagen specific Tregs with AAV8 suppresses autoimmune responses and arthritis in DBA1 mice through IL10 production. *Sci Rep*. **11**(1):18204. doi: 10.1038/s41598-021-97739-w.
20. Smith, C.M. et al. (2021). Phosphatidylserine binding directly regulates TIM-3 function. *Biochem J*. **478**(17):3331-3349. doi: 10.1042/BCJ20210425.
21. Kiss, E. et al. (2021). Artesunate restores the levels of inhibitory synapse proteins and reduces amyloid- $\beta$  and C-terminal fragments (CTFs) of the amyloid precursor protein in an AD-mouse model. *Mol Cell Neurosci*. doi: 10.1016/j.mcn.2021.103624.
22. Kiss, E. et al. (2021). Artemisinin-treatment in pre-symptomatic APP-PS1 mice increases gephyrin phosphorylation at Ser270: a modification regulating postsynaptic GABAAR density. *Biol Chem*. doi: 10.1515/hsz-2021-0153.
23. Jeong, J. et al. (2021). Combination of cell signaling molecules can facilitate MYOD1-mediated myogenic transdifferentiation of pig fibroblasts. *J Anim Sci Biotechnol*. **12**(1):64. doi: 10.1186/s40104-021-00583-1.
24. Martin-Jaular, L. et al. (2021). Unbiased proteomic profiling of host cell extracellular vesicle composition and dynamics upon HIV-1 infection. *EMBO J*. doi: 10.15252/embj.2020105492.
25. Xie, W. et al. (2021). Pterostilbene accelerates wound healing by modulating diabetes-induced estrogen receptor  $\beta$  suppression in hematopoietic stem cells. *Burns Trauma*. doi: 10.1093/burnst/tkaa045.
26. Sinha, A. et al. (2020). RAD51-Mediated DNA Homologous Recombination Is Independent of PTEN Mutational Status. *Cancers (Basel)*. **12**(11):E3178. doi: 10.3390/cancers12113178.
27. Sun, W. et al. (2020). snRNA-seq reveals a subpopulation of adipocytes that regulates thermogenesis. *Nature*. **587**(7832):98-102. doi: 10.1038/s41586-020-2856-x.
28. Yoshihara, E. et al. (2020). Immune-evasive human islet-like organoids ameliorate diabetes. *Nature*. doi: 10.1038/s41586-020-2631-z.
29. Baci, D. et al. (2020). Extracellular Vesicles from Skeletal Muscle Cells Efficiently Promote Myogenesis in Induced Pluripotent Stem Cells. *Cells*. **9**(6):E1527. doi: 10.3390/cells9061527.
30. Smolková, K. et al. (2020). SIRT3 and GCN5L regulation of NADP $^{+}$ - and NADPH-driven reactions of mitochondrial isocitrate dehydrogenase IDH2. *Sci Rep*. **10**(1):8677. doi: 10.1038/s41598-020-65351-z.

### **Warranty**

These products are warranted to perform as described in their labeling and in Cell Biolabs literature when used in accordance with their instructions. THERE ARE NO WARRANTIES THAT EXTEND BEYOND THIS EXPRESSED WARRANTY AND CELL BIOLABS DISCLAIMS ANY IMPLIED WARRANTY OF MERCHANTABILITY OR WARRANTY OF FITNESS FOR PARTICULAR PURPOSE. CELL BIOLABS's sole obligation and purchaser's exclusive remedy for breach of this warranty shall be, at the option of CELL BIOLABS, to repair or replace the products. In no event shall CELL BIOLABS be liable for any proximate, incidental or consequential damages in connection with the products.

***This product is for RESEARCH USE ONLY; not for use in diagnostic procedures.***

## **Contact Information**

Cell Biolabs, Inc.  
5628 Copley Drive  
San Diego, CA 92111  
Worldwide: +1 858-271-6500  
USA Toll-Free: 1-888-CBL-0505  
E-mail: [tech@cellbiolabs.com](mailto:tech@cellbiolabs.com)  
[www.cellbiolabs.com](http://www.cellbiolabs.com)

©2006-2026: Cell Biolabs, Inc. - All rights reserved. No part of these works may be reproduced in any form without permissions in writing.